

Expert on Flash Chromatography

Automated Flash Chromatography Systems



Smart Flash EPCLC AI-580S
Single Channel, 0-80ml/min

Intelligent & Eco-Friendly!

Fast and Efficient Separation

A target compound will always elute at around 4 column volume regardless of column size and flow rate, minimizing the run time while achieving a good separation.



Minimize Solvent Consumption.



Eco-Friendly!



Smart Flash EPCLC W-Prep 2XY
Dual Channel, 0-80ml/min
More Compact, Modular & Space Saving



Smart Flash EPCLC AI-700X
Single Channel, 5-150ml/min
5-200ml/min
15-400ml/min