


**Expert on Flash Chromatography**

# Automated Flash Chromatography Systems



**“New” AI-580**  
Single Channel, 0-80ml/min  
Easy Access to Solvent Bottles

## Intelligent & Eco-Friendly!

### Fast and Efficient Separation

A target compound will always elute at around 4 column volume regardless of column size and flow rate, minimizing the run time while achieving a good separation.



**Minimize Solvent Consumption.**



**Eco-Friendly!**



**W-Prep 2XY A Type**  
Dual Channel, 0-80ml/min  
More Compact, Modular & Space Saving



**W-Prep 2XY B Type**  
Dual Channel, 0-80ml/min  
Easy & Tight Column Connection



YAMAZEN CORPORATION