


“Universal Premium”, World’s Best Disposable Column




- Able to take double the sample as that of an ordinary column.
Universal 3L Premium (30 μ m, Silica, 135g)



| | |
|--------------------------|--------|
| Sample: | |
| Butyl-p-hydroxybenzoate | 675mg |
| Methyl-p-hydroxybenzoate | 675mg |
| Total. | 1350mg |


A high performance column free of peak fronting and peak tailing assures excellent sample purification and separation, even when a large amount of sample is loaded.

If a column produces fronting peaks and/or tailing peaks, then the amount of sample that a column can take will have to be radically reduced.



- Downsizing from L-size column to M-size column

Hi-Flash L (40 μ m, Silica, 30g)



Universal M (30 μ m, Silica, 16g)

