


## Further Enhancement of Separation Performance by using Inject Column Premium (40-63 $\mu$ m) and Universal Column Premium (25-40 $\mu$ m)



Inject Column L-size

Inject Column L-size (70 $\mu$ m) + Universal Column Premium L-size (25-40 $\mu$ m, 40g)




Sample: Methyl-p-hydroxybenzoate  
Nifedipine




Hexane: Ethyl acetate = 70:30

Universal Column Premium L-size



**A Drastic Enhancement of Separation Performance!**

Inject Column Premium L-size (40-63 $\mu$ m) + Universal Column Premium L-size (25-40 $\mu$ m, 40g)



Inject Column Premium L-size



Universal Column Premium L-size



**Downsizing effect!  
Reduction of sample running cost!**

Hi-Flash Column Premium 3L-size



Hi-Flash Column Premium 3L-size (25-40 $\mu$ m, 120g)

